Portables



Portavo

Portavo 904 pH

Portable, sturdy process analyzer for pH and ORP measurement.

Up to 5,000 values can be recorded using the integrated data logger. Using the USB port and Paraly SW 112 software, the logger data can be easily transferred to a PC for evaluation.

Custom pH Calibration Cal SOP

The Cal SOP calibration procedure allows pH sensors to be checked with up to 3 calibration points. A buffer is used as the verification buffer. The buffer set for each calibration point can be separately selected, thus also allowing their order to be determined.

Custom buffer solutions can be used, or choose from a list of commercially available buffer solutions, e.g., CaliMat, NIST, and DIN. A maximum permissible deviation (Delta pH) is entered for the verification buffer.

Security Package, Including User Management

Professional user management regulates access to the device and the sensor.

- Increased security for configuration, calibration, measurement data, and data logger settings.
- No unauthorized interventions during the operating cycle
- Up to 4 user profiles can be set
- Different access rights can be established

Depending on the user's experience, the role profile can optionally be defined for configuration of the device and sensor or for sensor calibration. This clearly minimizes the risk of inadvertently changing settings.

Greater Reliability During Operation

Memosens sensors can be assigned directly to the Portavo using the data stored in the sensor, such as

Sensor type TAG Group

Unambiguous assignment of the sensor to the device reduces the potential for errors. This ensures that only the right sensors are used for the selected measuring point.




Knick >

pH/ORP Measurement



Facts and Features


- Memosens sensors or analog sensors for pH or ORP measurement can be used with one device
- Sensor quiver protects the sensor from drying out and damage
- Sturdy housing with IP66/67 protection, also suitable for outdoor use
- Li-ion rechargeable battery
 USB chargeable
- Data logger with 5,000 values
- Micro USB port and Paraly
 SW 112 operating software
- Custom pH calibration Cal SOP
- User management for access control
- Sensor verification for clear assignment of the sensor to the device via sensor type, TAG, or group
- Temperature detector adjustment in the Memosens sensor (offset correction)





A122604

061



meas	cal	set



Portables

Portavo

Specifications					
pH/mV input (analog)	pH socket DIN 19 262 (13/4	mm)			
	pH measuring range				
	Decimal places*)	2 or 3			
	Input resistance	1 x 10 ¹² Ω (0 35 °C / 32 86 °F)			
	Input current	1 x 10 ⁻¹² A (at RT, doubles every 10 K)			
	Measuring cycle Measurement error ^{1,2,3)}	Approx. 1 s			
		< 0.01 pH –1300 1300 mV	TC < 0.001 pH/K		
	mV measuring range Measuring cycle	Approx. 1 s			
	Measurement error ^{1,2,3)}	< 0.1 % of meas.	TC < 0.03 mV/K		
		value + 0.3 mV			
Temperature input	2 x Ø 4 mm for integrated or separate temperature detector				
	Measuring ranges	NTC 30 kΩ	–20 120 °C / –4 248 °F		
	5 5	Pt1000	–40 250 °C / –40 482 °F		
	Measuring cycle	Approx. 1 s			
	Measurement error ^{1,2,3)}	< 0.2 K (Tamb = 23 °C / 73.4 °F); TC $<$ 25 ppm/K			
Memosens pH input, ISFET	M8 socket, 4-pin, for Memo	M8 socket, 4-pin, for Memosens laboratory cable			
	Display ranges ⁴⁾	рН	-2.000 16.000		
		mV	–2000 2000 mV		
		Temperature	–50 250 °C / –58 482 °F		
Sensor adjustment*)	pH calibration				
Operating modes ^{*)}	Calimatic	Calibration with automatic buffer recognition			
	Manual	Manual calibration with entry of individual buffer values			
	ISFET zero	Calibration of ISFET sensors			
	Data entry	Data entry of zero and slope			
	Cal-SOP (TAN option)	Software option SW-P001: Defining the pH buffers and the sequence of the calibration steps; defining the delta deviation for the verification buffer			
	Temperature calibration (TAN option)	Software option SW-P002 for temperature detec- tor adjustment in the Memosens sensor (offset correction)			
Calimatic buffer sets*)	–01– Mettler-Toledo	2.00/4.01/7.00/9.2	1		
	–02– Knick CaliMat	2.00/4.00/7.00/9.00	0/12.00		
	–03– Ciba (94)	2.06/4.00/7.00/10.0	2.06/4.00/7.00/10.00		
	–04– NIST Technical	1.68/4.00/7.00/10.01/12.46			
	–05– NIST Standard	1.679/4.006/6.865/9.180			
	-06- HACH		4.01/7.00/10.01/12.00		
	-07- WTW techn. buffers	2.00/4.01/7.00/10.00			
	–08– Hamilton		2.00/4.01/7.00/10.01/12.00		
	-09- Reagecon	2.00/4.00/7.00/9.00/12.00			
	–10– DIN 19267 –U1– (User)	1.09/4.65/6.79/9.23/12.75			
Permissible calibration range		Loadable via Paraly SW 112			
	Zero point Slope	6 8 pH			
	ISFET	Approx. 74 104 % -750 +750 mV Operating point (asymmetry)			
	Calibration timer ^{*)}				
		Interval 1 99 days, can be deactivated Provides information on the condition of the sense			
	Evaluation of				
	EVALUATION OF	Zero point/slope, response time, cal. interval			



pH/ORP Measurement

Specifications

Memosens ORP input	M8 socket, 4-pin, for Memosens laboratory cable				
	Display ranges ⁴⁾	mV Temperature	–2000 2000 mV –50 250 °C / −58 482 °F		
Sensor adjustment*)	ORP calibration (zero offset) Permissible calibration range	ΔmV (offset)	–700 700 mV		
	Temperature calibration (TAN option)	Software option SW-P002 for temperature detec- tor adjustment in the Memosens sensor (offset correction)			
Connections	2 x socket Ø 4 mm for separate temperature probe 1 x M8 socket, 4-pin, for Memosens laboratory cable 1 x micro USB-B for data transmission to PC 1 x pH socket in acc. with DIN 19262				
Display	LCD STN 7-segment display w Status indicators Notices	rith 3 lines and icons For battery condition, logger Hourglass			
Keypad		[on/off], [cal], [meas], [set], [▲], [▼], [STO], [RCL], [clock]			
Data logger	Space for 5,000 entries Recording	Manual, interval- or event-controlled			
MemoLog calibration data logger (Memosens only)	Can save up to 100 Memosen Recording	s calibration records Directly readable via MemoSuite or Paraly SW 112 (USB)			
	Can be shown on the display	Manufacturer, sensor type, serial no., zero point, slope, calibration date			
Communication	USB 2.0				
	Profile Usage	HID, driverless installation Data transfer and configuration via the Paraly SW 112 software			
Diagnostic functions	Sensor data (Memosens only)) Manufacturer, sensor type, serial number, operating time			
	Calibration data Device self-test Device data	Calibration date; zero point, slope Automatic memory test (FLASH, EEPROM, RAM) Device type, software version, hardware version			
Data retention	Parameter, calibration data > 10 years				
EMC	EN 61326-1 (General requiren Emitted interference	,	-0		
	Immunity to interference EN 61326-2-3	Class B (residential) Industrial applications			
RoHS conformity	According to Directive 2011/65/EU				
Power supply	4 x AA (Mignon) alkaline batteries 4 x NiMH rechargeable batteries or 1 x Li-ion rechargeable battery (rechargeable via USB) Operating time Approx. 1000 h (alkaline)				
Rated operating conditions	Ambient temperature Transport / storage temp. Relative humidity	-10 55 °C / 14 -25 70 °C / −13 0 95 %, brief co			
Housing	Material Ingress protection Dimensions Weight	PA12 GF30 + TPE IP66/67 with pressure compensation Approx. 132 x 156 x 30 mm / 5.2 x 6.14 x 1.18 inche Approx. 500 g / 1.10 lbs			
*) User-defined	2) + 1 digit				

At rated operating conditions
 Plus sensor error

2) ± 1 digit4) Ranges dependent on Memosens sensor