




# neodisher<sup>®</sup> LM 2



## Detergent for the manual reprocessing of thermostable and thermolabile instruments and laboratory glassware



Liquid concentrate

### Fields of application:

- Manual cleaning of laboratory glassware and surgical instruments in immersion and ultrasonic baths in hospitals, dental practices, medical, industrial and biological laboratories
- Manual cleaning of utensils and animal cages
- Cleaning of pipettes in special washers

### Performance spectrum:

- Removes organic residues, e.g. encrusted blood and protein residues as well as feed residues and droppings from animal cages and inorganic residues
- Supports the removal of dental cement
- Easily rinsed off
- Does not produce any interfering foam
- Also suitable for manual pre-cleaning due to the low-foam formulation the items to be cleaned can subsequently be washed in the machine without without prior rinsing
- Suitable for commonly used laboratory utensil materials, e.g. glass, ceramic, stainless steel and plastic
- Objects made of light metals should be tested first
- Usable at all degrees of water hardness

### Special properties:

- Contains low-foaming surfactants and various cleaning enhancers
- Especially suitable for the ultrasonic cleaning of surgical and dental instruments, particularly micro instruments
- If laboratory glassware is autoclaved before machine cleaning, a wet autoclaving with the addition of neodisher LM 2 (30 ml/l) is recommended. Afterwards residues of soiling can be more easily removed during automated cleaning

### Application and dosage:

Manual cleaning of laboratory glassware and surgical instruments:  
5 - 20 ml/l in immersion and ultrasonic baths  
Contact time: 10 - 30 min each

Manual cleaning of utensils and animal cages:  
5 – 30 ml/l  
Contact time: 10 – 30 min

Automated cleaning of pipettes:  
10 - 30 ml/l in pipette washers,  
Contact time: 10 - 30 min

Wet autoclaving of laboratory glassware:  
Add 30 ml/l neodisher LM 2. Dirt residues do not burn in and can be removed more easily during subsequent machine cleaning.

The neodisher LM 2 solution has to be rinsed off completely with water (preferably deionised water).

### Notes on application:

- Do not mix with other products
- For professional use only
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664

### Technical data:

pH value	10.6 - 11.3 (5 - 30 ml/l, determined in deionised water, 20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)
Density	approx. 1.2 g/cm <sup>3</sup> (20 °C)

### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:  
< 5 % amphoteric surfactants, anionic surfactants  
15 - 30 % phosphates




# neodisher<sup>®</sup> LM 2

## CE-mark:

neodisher LM<sub>2</sub> complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

The product is sensitive to frost. Always store at a temperature between 0 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4403/3-1

Date of issue: 02/2021