




Scientific Equipment

Tel. 051 501153 - Fax 051 6336182

[www.favs.it](http://www.favs.it) - [info@favs.it](mailto:info@favs.it)

# IKA

designed for scientists



## R 1331 Anchor stirrer

/// Data Sheet

Tangential flow, high shearing rate at edges, minimum deposits on the vessel wall. Used at low speeds. Polymer reactions, even distribution of high mineral contents in liquids. The ideal stirrer\* for medium to highly viscous fluids.

Stirrer Ø: 90 mm

Shaft Ø: 8 mm

Shaft length: 350 mm

[www.ika.com](http://www.ika.com)

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

### Technical Data

|                                 |                    |
|---------------------------------|--------------------|
| Stirrer diameter [mm]           | 90                 |
| Shaft diameter [mm]             | 8                  |
| Shaft length [mm]               | 350                |
| Speed max. [rpm]                | 1000               |
| Material in contact with medium | AISI 316L / 1.4571 |
| Weight [kg]                     | 0.2025             |

