



Model:
Heating stage

Description:
Heating stage for temperature control of specimens

FAVS
Scientific Equipment

Tel. 051501153

www.favs.it • info@favs.it

Stage	The heating elements are mounted directly inside the stage of the microscope [refer to the datasheet of your microscope model for further information].
Power supply	External Power Supply. Input 100/240 Vac, 50/60 Hz; Output 6 Vdc 5 A. Max. required power: 25 W.
Temperature specifications	<ul style="list-style-type: none"> • Temperature controlling range: from 20° C up to 50° C. • Temperature setting resolution: 0.5° C. • Temperature measurement resolution: 0.1° C. • Average error between actual stage temperature and displayed temperature: +/- 1° C.
External controller	<ul style="list-style-type: none"> • 16 characters, back-illuminated LCD display. • Rotating knob for temperature adjusting. • ON/OFF switch.