



















# Cuvette

<b>Vetro ottico 320 – 1400 nm</b>						
<b>Cammino ottico (mm)</b>	1	10	20	30	50	100
<b>Dimensioni esterne WxLxH (mm)</b>	12.5x3.5x45	12.5x12.5x45	12.5x22.5x45	12.5x32.5x45	12.5x52.5x45	12.5x102.5x45
<b>Volume (ml)</b>	0.4	3.5	7.0	10.5	17.5	35.0
<b>Facce ottiche</b>	2	2	2	2	2	2
<b>Tappo in Teflon</b>	✓	✓	✓	✓	✓	✓
<b>Numero di cuvette</b>	2	4	4	4	4	4
<b>Codice</b>	<b>11000382</b>	<b>11000252</b>	<b>11000262</b>	<b>11000272</b>	<b>11000282</b>	<b>11000292</b>

<b>Quarzo 190 – 1400 nm</b>						
<b>Cammino ottico (mm)</b>	1	10	20	30	50	100
<b>Dimensioni esterne WxLxH (mm)</b>	12.5x3.5x45	12.5x12.5x45	12.5x22.5x45	12.5x32.5x45	12.5x52.5x45	12.5x102.5x45
<b>Volume (ml)</b>	0.4	3.5	7.0	10.5	17.5	35.0
<b>Facce ottiche</b>	2	2	2	2	2	2
<b>Tappo in Teflon</b>	✓	✓	✓	✓	✓	✓
<b>Numero di cuvette</b>	2	2	2	2	2	2
<b>Codice</b>	<b>11000522</b>	<b>11000302</b>	<b>11000312</b>	<b>11000322</b>	<b>11000332</b>	<b>11000342</b>

<b>Celle speciali</b>						
<b>Cammino ottico (mm)</b>	10	10	10	10	10	10
<b>Dimensioni esterne WxLxH (mm)</b>	12.5x12.5x45	12.5x12.5x45	12.5x12.5x45	12.5x12.5x45	12.5x12.5x45	12.5x12.5x45
<b>Altezza del raggio (mm)</b>	15	15	15	15	15	15
<b>Volume (µl)</b>	100	200	500	1000 (1ml)	700	700
<b>Materiale</b>	Quarzo	Quarzo	Quarzo	Quarzo	Vetro ottico	Quarzo
<b>Range (nm)</b>	190 -1400	190 -1400	190 -1400	190 -1400	320 -1400	190 -1400
<b>Cella a flusso</b>	–	–	–	–	✓	✓
<b>Facce ottiche</b>	2	2	2	2	2	2
<b>Tappo in Teflon</b>	✓	✓	✓	✓	–	–
<b>Numero di cuvette</b>	1	1	1	1	1	1
<b>Codice</b>	<b>11000352</b>	<b>11000362</b>	<b>11000372</b>	<b>11000392</b>	<b>11000492</b>	<b>11000502</b>