



# MST

## Magnetic Stirrers

---

Small and efficient magnetic stirrer



# Resistant and Easy-to-use

With its white surface the MST facilitates colour change detection point making it suitable for many applications including microbiology and biochemistry, microtitration stirrings and slow/long stirring.

- Speed regulation up to 1100 rpm
- Accepts up to 5 liter flasks
- Compact design
- Technopolymer casing resists chemicals
- Excellent speed control, even at low speeds
- Ideal for microtitration, BOD, microbiology and biochemistry



## INSTRUMENT - CODE

MST 100-240 V / 50-60 Hz F203A0440

## TECHNICAL DATA

|                             |                                    |
|-----------------------------|------------------------------------|
| CONSTRUCTION MATERIAL       | Technopolymer                      |
| STIRRING VOLUME             | Up to 5 L                          |
| ELECTRONIC SPEED REGULATION | Up to 1100 rpm                     |
| STIRRING SYSTEM             | High power magnet type PCM         |
| DIMENSION (WXHXD)           | 120x50x145 mm<br>(4.7x1.97x5.7 in) |
| POWER                       | 0.6 W                              |
| WEIGHT                      | 0.4 Kg (0.9 lb)                    |



We reserve the right to make technical alterations  
We do not assume liability for errors in printing, typing or transmission

Rev 1 04.2021

**favs**  
Scientific Equipment  
Tel. 051 501153 - Fax 051 6336182  
[www.favs.it](http://www.favs.it) - [info@favs.it](mailto:info@favs.it)

**DESIGNED AND MANUFACTURED IN ITALY**

**VELP SCIENTIFICA**

|  |   |  |
|--|---|--|
| <b>ITALY - HQ</b><br>Via Stazione 16<br>20865 Usmate (MB) Italy<br>Tel. +39 039 628811<br><a href="mailto:velpitalia@velp.com">velpitalia@velp.com</a> | <b>USA</b><br>155 Keyland Court, Bohemia<br>NY 11716 - U.S.<br>Tel. +1 631 573 6002<br><a href="mailto:velpusa@velp.com">velpusa@velp.com</a> | <b>CHINA</b><br>Xinlong Rd Building 28, Lane 1333<br>Shanghai city - China<br>Tel. +86 180 1755 7329<br><a href="mailto:velpchina@velp.com">velpchina@velp.com</a> |
| <b>INDIA</b><br><a href="mailto:velpindia@velp.com">velpindia@velp.com</a>   |   |  |