


### ▶ Horizontal electrophoresis cell

- 1 | Bucket made of high quality acrylic material and in one piece to prevent liquid leakage.
- 2 | Transparent cover that provides a correct seal to avoid buffer evaporation and electrical leakage; with fixed power cables with banana type connector.
- 3 | Easy to replace platinum electrodes.
- 4 | The volume of buffer required guarantees not only the refrigeration but also keeps the pH values stable throughout the process.
- 5 | Auto-off function when lid is removed.

### ▶ INCLUDED ACCESSORIES

- 1 | Mold for the preparation of gels
- 2 | Dark-band gel tray for easy sample loading
- 3 | Set of combs of different thickness and number of samples.



Code	ZFD013
Gel size	78x100 mm
Number of samples	8/15
Combs size	1 y 2 mm
Buffer (approx.)	300 mL aprox.
Input voltage	150 V
Dimensions (LxWxH)	365x105x55 mm

## HORIZONTAL ELECTROPHORESIS CELL [Z] OTHERS